

## Product datasheet for **RC206740L2V**

### **HID1 (NM\_030630) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	HID1 (NM_030630) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HID1
Synonyms:	17orf28; C17orf28; DMC1; HID-1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_030630
ORF Size:	2364 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206740).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_030630.1</a> , <a href="#">NP_085133.1</a>
RefSeq Size:	3369 bp
RefSeq ORF:	2367 bp
Locus ID:	283987
UniProt ID:	<a href="#">Q8IV36</a>
Cytogenetics:	17q25.1
MW:	88.7 kDa



[View online »](#)

**Gene Summary:** May play an important role in the development of cancers in a broad range of tissues.  
[UniProtKB/Swiss-Prot Function]