

Product datasheet for **RC206749A1V**

Human C1orf83 (TCEANC2) (NM_153035) AAV Particle

Product data:

Product Type: AAV Particles
 Product Name: Human C1orf83 (TCEANC2) (NM_153035) AAV Particle
 Tag: Myc-DDK
 Symbol: C1orf83
 Synonyms: C1orf83
 Mammalian Cell Selection: None
 Vector: pAAV-AC-Myc-DDK (PS100089)
 ORF Nucleotide Sequence: >RC206749 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGATAAATTCGTCAATCGAACGCCTAGAATCCAGAATAGCCCTCAGAAGAAAGATTCTGGAGGAAAGG
 TTTACAAGCAGGCCACGATTGAATCTCTGAAGAGAGTTGTGGTTGTAGAAGACATAAAAAAGATGAAAAAC
 TATGCTGGAGCTTCCTGATCAACCAAGAGAAATCTTGTGAAGCCTTACAAGAATTAAGAAGAAAAATA
 CCCTCCAGGGAAGTGTTAAAAATCAACAAGGATAGGTCACACTGTGAACAAGATGCGTAAACACTCAGATT
 CAGAAGTGGCTTCTCTTGCCAGAGAAGTTTACACTGAGTGGAAAACTTTCACTGAAAAACATTCAAATAG
 ACCTTCTATTGAAGTTAGAAGTGATCCCAAAACCGAGTCGTTGAGGAAAAATGCTCAGAAATTACTCTCA
 GAAGCCTTGAATTAAGATGGATCACCTACTGTTGAAAATATTGAACGGGAAACGTTTCATCTCTGCT
 CCCGCCTCATTAATGGGCCGTACCGGCGGACGGTGAGAGCCCTGGTCTTCACATTAAGCACCGAGCTGA
 AATCCGGGCTCAGGTGAAGAGCGGCTCGCTGCCAGTCGGCACGTTTGTACAGACCCACAAAAAG

ACGCGTACGCGGCGGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC206749 protein sequence
 Red=Cloning site Green=Tags(s)

MDKFVIRTPRIQNSPQKKDSGGKVYKQATIESLKRVVVEDIKRWKTMLELPDQTKENLVEALQELKKKI
 PSREVLKSTRIGHTVNKMRKHSDEVASLAREVYTEWKTFTKHSNRPSIEVRSDPKTESLRKNAQKLLS
 EALELKMDHLLVENIERETFHLCSRLINGPYRRTVRALVFTLKHRAEIRAQVKSGLPVGTFVQTHKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_153035
ORF Size:	624 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_153035.1</u>
RefSeq Size:	4097 bp
RefSeq ORF:	627 bp
Locus ID:	127428
UniProt ID:	<u>Q96MN5</u>
Cytogenetics:	1p32.3
MW:	24.2 kDa