

Product datasheet for **RC206753A1V**

Human NR2F2 (NM_021005) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human NR2F2 (NM_021005) AAV Particle
Tag:	Myc-DDK
Symbol:	NR2F2
Synonyms:	ARP-1; ARP1; CHTD4; COUPTF2; COUPTFB; COUPTFII; NF-E3; SRXX5; SVP40; TFCOUP2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206753 representing NM_021005 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCAATGGTAGTCAGCACGTGGCGGACCCCCAGGACGAGGTGCCCGGCTCACAGGGCAGCCAGGCCT
CGCAGGCGCGCCCGTGGCCGGCCCGCCGCGCCCGCACACGCCACAGACGCCCGGCCAAGGGGG
CCCAGCCAGCAGCCAGCCAGCCAGCGCGCCGGTGGCCAGGGCGGCCCTGGCGGCCCGGGTAGCGACAAG
CAGCAGCAGCAGCAACACATCGAGTGCCTGGTGTGCGGAGACAAGTCGAGCGGCAAGCACTACGCCAGT
TCACGTGCGAGGGCTGCAAGAGCTTCTTCAAGCGCAGCGTGCGGAGGAACCTGAGCTACACGTGCCGCGC
CAACCGGAACGTGCCATCGACCAGCACCATCGCAACCAAGTCCAGTACTGCCGCTCAAAAAGTGCCTC
AAAGTGGGATGAGACGGAAGCGGTGCAGAGGGGAGGATGCCGCGACCCAGCCAGCCAGCCAGGGCAGT
TCGCGCTGACCAACGGGGATCCCCTCAACTGCCACTCGTACCTGTCCGGATATATTTCCCTGCTGTTGCG
CGCGGAGCCCTATCCACGTGCGCTTCGGCAGCCAATGCATGCAGCCCAACAACATCATGGGTATCGAG
AACATTTGCGAACTGGCCGCGAGGATGCTCTTCAGCGCCGTCGAGTGGGCCGGAACATCCCCTTCTTCC
CCGACCTGCAGATCACGACAGGTGGCCCTGCTTCGCTCACCTGGAGCGAGCTGTTTGTGTTGAATGC
GGCGCAGTGTCCATGCCCTCCACGTGCCCCGCTCCTGGCCGCCCGGCCCTGCATGCTTCGCCCATG
TCCGCCGACCGGGTGGTGCCTTTATGGACCACATACGGATCTTCCAAGAGCAAGTGGAGAAGCTCAAGG
CGCTGCAGTTGACTCAGCCGAGTACAGCTGCCTCAAGGCCATAGTCCTGTTACCTCAGATGCCTGTGG
TCTCTCTGATGTAGCCATGTGGAAGCTTGCAGGAAAAGTCTCAGTGTGCTTTGGAAGAATACGTTAGG
AGCCAGTACCCCAACAGCCGACGAGATTGGAAGGCTTTTGCTTCGCTCCCTCCCTCCGACCGTCT
CCTCCTCAGTCATAGAGCAATTGTTTTTCGTCCGTTTGGTAGGTAACCCCATCGAAACCTCATCCG
GGATATGTTACTGTCCGCAGCAGTTTTAACTGGCCGTATATGGCAATTCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206753 representing NM_021005 Red=Cloning site Green=Tags(s) MAMVVSTWRDPQDEVPGSQSQASQAPPVPGPPPGAPHTPQTPGQGGPASTPAQTAAAGGQGGPGGPGSDK QQQQHIIECVVCGDKSSGKHYGQFTCEGCKSFFKRSVRRNL SYTCRANRNCPIQHHRNQCQYCRLLKKCL KVGMRREAVQGRMPPTQPTHGQFALTNGDPLNCHSYLSGYISLLLRAEPYPTSRFGSQCMQPNNIMGIE NICELAAARMLFSAVEWARNIPFFPDQLITDQVALLRLTWSELFVLNAAQCSMPLHVAPLLAAAGLHASPM SADRVVAFMDHIRIFQEYQVEKLKALHVDSA EYSCLKAIVLFTSDACGLSDVAHVESLQEKSQCALEEYVR SQYPNQPTRFGKLLLRPSLRTVSSSVIEQLFFVRLVGKTP IETLIRDMLLSGSSFNWPYMAIQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_021005
ORF Size:	1242 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_021005.2
RefSeq Size:	3049 bp
RefSeq ORF:	1245 bp
Locus ID:	7026
UniProt ID:	P24468
Cytogenetics:	15q26.2
MW:	45.4 kDa