

Product datasheet for **RC206762L3V**

GABARAPL1 (NM_031412) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GABARAPL1 (NM_031412) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GABARAPL1
Synonyms:	APG8-LIKE; APG8L; ATG8; ATG8B; ATG8L; GEC1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031412
ORF Size:	351 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206762).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031412.2
RefSeq Size:	1885 bp
RefSeq ORF:	354 bp
Locus ID:	23710
UniProt ID:	Q9H0R8
Cytogenetics:	12p13.2
Domains:	MAP1_LC3
Protein Pathways:	Regulation of autophagy



[View online »](#)

MW: 14 kDa

Gene Summary: Ubiquitin-like modifier that increases cell-surface expression of kappa-type opioid receptor through facilitating anterograde intracellular trafficking of the receptor. Involved in formation of autophagosomal vacuoles. Whereas LC3s are involved in elongation of the phagophore membrane, the GABARAP/GATE-16 subfamily is essential for a later stage in autophagosome maturation.[UniProtKB/Swiss-Prot Function]