

Product datasheet for **RC206773L3V**

Membrin (GOSR2) (NM_001012511) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Membrin (GOSR2) (NM_001012511) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Membrin
Synonyms:	Bos1; EPM6; GS27
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001012511
ORF Size:	585 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206773).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001012511.1
RefSeq Size:	2220 bp
RefSeq ORF:	588 bp
Locus ID:	9570
UniProt ID:	O14653
Cytogenetics:	17q21.32
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	SNARE interactions in vesicular transport



[View online »](#)

MW: 22.7 kDa

Gene Summary: This gene encodes a trafficking membrane protein which transports proteins among the medial- and trans-Golgi compartments. Due to its chromosomal location and trafficking function, this gene may be involved in familial essential hypertension. [provided by RefSeq, Mar 2016]