

## Product datasheet for **RC206778L2V**

### **CYTL1 (NM\_018659) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	CYTL1 (NM_018659) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CYTL1
Synonyms:	C4orf4; C17
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_018659
ORF Size:	408 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206778).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018659.2</a>
RefSeq Size:	1019 bp
RefSeq ORF:	411 bp
Locus ID:	54360
UniProt ID:	<a href="#">Q9NRR1</a>
Cytogenetics:	4p16.2
Protein Families:	Druggable Genome, Secreted Protein
MW:	15.6 kDa



[View online »](#)

**Gene Summary:**

C17 is a cytokine-like protein specifically expressed in bone marrow and cord blood mononuclear cells that bear the CD34 (MIM 142230) surface marker (Liu et al., 2000 [PubMed 10857752]).[supplied by OMIM, Mar 2008]