

Product datasheet for **RC206785L4V**

PRR15 (NM_175887) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PRR15 (NM_175887) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PRR15
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_175887
ORF Size:	387 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206785).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_175887.2
RefSeq Size:	1717 bp
RefSeq ORF:	390 bp
Locus ID:	222171
UniProt ID:	Q8IV56
Cytogenetics:	7p14.3
MW:	13.7 kDa
Gene Summary:	May have a role in proliferation and/or differentiation.[UniProtKB/Swiss-Prot Function]


[View online »](#)