

Product datasheet for **RC206792A1V**

Human TAF4 (NM_182522) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human TAF4 (NM_182522) AAV Particle
Tag: Myc-DDK
Symbol: TAF4
Synonyms: FAM19A4; TAF4-4
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC206792 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAGGTCCCCAAGGATGAGAGTCTGTGCTAAGTCAGTGTGCTGTCGCACTGGCTCTTTCTAGCCTACG
 TGTTAATGGTGTGCTGTAAGCTGATGTCCGCTCAAGCCAGCACCTCCGGGGACATGCAGGTCACCACCA
 AATCAAGCAAGGGACCTGTGAGGTGGTCGCGTGCACAGGTGCTGCAATAAGAACCGCATAGAAGAGCGG
 TCACAAACGGTCAAGTGCTCTTGCTTCCCGGGACAGGTGGCGGGCACAACCTCGGGCTCAACCTTCTTG
 TTGAAGCTTCCATTGTGATTAGAAATGGTGGTGTACATGAATCCGTGTTGGAAGGAGAGGATTGTAA
 AGTGCTGCCAGATTACTCAGTTGGTCTGTAGCAGTGGCAATAAAGTCAAACTACGAAGGTAACGCGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC206792 protein sequence
 Red=Cloning site Green=Tags(s)

MRSPMRVCAKSVLLSHWLFLAYVLMVCCKLMSASSQHLRGHAGHHQIKQGTCEVVAVHRCCKNKNRIEER
 SQTVKCSCFPQVAGTTTRAQPSVEASIVIQWCHMNPCLGEDCKVLPDYSWSSSGNKVKTTKVTR

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_182522



[View online »](#)

ORF Size:	420 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_182522.3</u>
RefSeq Size:	2294 bp
RefSeq ORF:	423 bp
Locus ID:	151647
UniProt ID:	<u>Q96LR4</u>
Cytogenetics:	3p14.1
MW:	15.7 kDa