

## Product datasheet for **RC206798L4V**

### DEDD2 (NM\_133328) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	DEDD2 (NM_133328) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DEDD2
Synonyms:	FLAME-3; FLAME3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133328
ORF Size:	978 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206798).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_133328.2</a> , <a href="#">NP_579874.1</a>
RefSeq Size:	2025 bp
RefSeq ORF:	981 bp
Locus ID:	162989
UniProt ID:	<a href="#">Q8WXF8</a>
Cytogenetics:	19q13.2
Domains:	DED
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**MW:** 36.2 kDa

**Gene Summary:** This gene encodes a nuclear-localized protein containing a death effector domain (DED). The encoded protein may regulate the trafficking of caspases and other proteins into the nucleus during death receptor-induced apoptosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2012]