

Product datasheet for **RC206800**

C11orf45 (NM_145013) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C11orf45 (NM_145013) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C11orf45
Synonyms: FLJ43646; MGC35558
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC206800 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTTGACGGCTGGTCTCAGTGCACACCTGAGTAGCACGACCTCTCCGCCCTGGACGCACGCTGCCA
 TCAGCTGGGAGCTGGACAACGTGCTGATGCCTAGTCCCAGAATCTGGCCCAGGTGACTCCAACAGGCAG
 GTCTGCCTCTGTCAGGAGTGAGGGTAACACCTCCTCACTCTGGAATTTCTCAGCTGGCAGGATGTGCAT
 GCCATAGTAACCAACCTGTGAGTCTGTGCTGAGCTCTGCCGTCTACACCCACGGCTGTGGCTGTGTGA
 GGTCTGCCACAAACATTACCTGTCAGTCTCAGGACAACAAAGGCAGGCCGCCCGCAGGAAGAGGAGAA
 CTCAATCTGCAAGGCCCATGATAGTAGAGGGCCGCTGGGCTACCCCTCAGTGCCCATCAGCCTGTT
 TCCGGTGGTCTAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC206800 protein sequence
 Red=Cloning site Green=Tags(s)

MLTRLVLSAHLSTTSPWTHAAISWELDNVLMSPRIWPQVPTGRSASVRSEGNTSSLWNFSAGQDVH
 AIVTRTCEVLSAVYTHGCGVRSATNITCQSSGQQRQAARQEEENSICKAHSREGRLGYPLSAHQPG
 SGGPN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms: https://cdn.origene.com/chromatograms/mk6015_a07.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_145013

ORF Size: 435 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_145013.2](#)


RefSeq Size: 3648 bp

RefSeq ORF: 438 bp

Locus ID: 219833
UniProt ID: [Q8TAV5](#)
Cytogenetics: 11q24.3
Protein Families: Secreted Protein
MW: 15.6 kDa

Product images:


Circular map for RC206800



Western blot validation of overexpression lysate (Cat# [LY403415]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC206800 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).