

Product datasheet for **RC206817**

FAM104A (NM_032837) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: FAM104A (NM_032837) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: FAM104A
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC206817 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCGAGCGCCTCCGTCCCAGGAAAAGGAGAAGGAATGGCAACGAAGAAGACAACCATCTTCCCCCCC
AGACAAAAGAAGTAGCAGAAACCCTGTCTTTCAGGATTCCTGGGACACAGAGTCTTCAGGCAGTGACAG
TGGTGGGAGCAGCAGCAGCAGCAGCAGCATCAACAGCCGGACAGGGCCAGCGGGCCGAAGGCAGC
TTGAGCCAGACCATGGCCGGATCCAGCCCTAACACGCCTCAGCCCGTGCCGAGCAGTCCGCGCTGTGCC
AAGGCCTCTACTTCCACATCAACCAGACCCTGAGGGAGGCCACTTCCACAGCCTACAGCACCGAGGGCG
GCCTCTGACA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC206817 protein sequence
Red=Cloning site Green=Tags(s)
MSERLRPRKRRRNGEEDNHLPPQTKRSSRNPFVQDSWDTESGSDSGGSSSSSSINSPPDRASGPEGS
LSQTMAGSSPNTQPVPPEQSALCQGLYFHINQTLREAHFHSLQHRGRPLT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6433_e12.zip



Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_032837

ORF Size: 360 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_032837.1](#), [NP_116226.1](#)


RefSeq Size: 2842 bp

RefSeq ORF: 561 bp


Locus ID: 84923

UniProt ID: Q969W3
 Cytogenetics: 17q25.1
 MW: 13.2 kDa

Product images:



Circular map for RC206817



Western blot validation of overexpression lysate (Cat# [LY409820]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC206817 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).