

## Product datasheet for RC206820L4

## ME1 (NM\_002395) Human Tagged Lenti ORF Clone

## **Product data:**

**Product Type:** Expression Plasmids

Product Name: ME1 (NM\_002395) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: ME1

Synonyms: HUMNDME; MES

Mammalian Cell Puromycin

Selection:

**Vector:** pLenti-C-mGFP-P2A-Puro (PS100093)

Sgfl-Mlul

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC206820).

Sequence:

equence.

Restriction Sites: Cloning Scheme:





<sup>\*</sup> The last codon before the Stop codon of the ORF

**ACCN:** NM\_002395

ORF Size: 1716 bp



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



**OTI Disclaimer:** 

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>

**OTI Annotation:** 

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 

- 1. Centrifuge at 5,000xg for 5min.
- 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
- 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** <u>NM 002395.3</u>

 RefSeq Size:
 3519 bp

 RefSeq ORF:
 1719 bp

 Locus ID:
 4199

 UniProt ID:
 P48163

 Cytogenetics:
 6q14.2

Domains: malic

**Protein Pathways:** Metabolic pathways, PPAR signaling pathway, Pyruvate metabolism

MW: 64.1 kDa

**Gene Summary:** This gene encodes a cytosolic, NADP-dependent enzyme that generates NADPH for fatty acid


biosynthesis. The activity of this enzyme, the reversible oxidative decarboxylation of malate, links the glycolytic and citric acid cycles. The regulation of expression for this gene is complex. Increased expression can result from elevated levels of thyroid hormones or by higher

icleased expression carries are from elevated levels of trigroup from ones of by fr


proportions of carbohydrates in the diet. [provided by RefSeq, Jul 2008]



## **Product images:**



Circular map for RC206820L4



Double digestion of RC206820L4 using Sgfl and Mlul  $\,$