

Product datasheet for **RC206831A1V**

Human SEPTIN5 (NM_002688) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SEPTIN5 (NM_002688) AAV Particle
Tag:	Myc-DDK
Symbol:	SEPTIN5
Synonyms:	CDCREL; CDCREL-1; CDCREL1; H5; HCDCREL-1; PNUTL1; SEPT5
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206831 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGCACAGGCCTGCGGTACAAGAGCAAGCTGGCGACCCAGAGGACAAGCAGGACATTGACAAGCAGT
ACGTGGGCTTCGCCCACTGCCCAACCAGGTGCACCGCAAGTCGGTGAAGAAAGGCTTTGACTTCACACT
CATGGTGGCTGGTGAGTCAGGCCTGGGGAAGTCCACACTGGTCCACAGCCTCTTCCTGACAGACTTGATC
AAGGACCGGAAGCTGCTCAGTGCTGAGGAGCGCATCAGCCAGACGGTAGAGATTCTAAACACACCGGTGG
ACATTGAGGAGAAGGGAGTCAAGCTGAAGCTCACCATCGTGGACACGCCGGGATTCCGGGACGCTGTCAA
CAACACCGAGTGCTGGAAGCCCATCACCGACTATGTGGACAGCAGTTTGAGCAGTACTCCCGTGATGAG
AGCGGCCTCAACCGAAAGAACATCCAAGACAACCGAGTGCACTGCTGCCTATACTTCATCTCCCCCTTCG
GGCATGGGCTGCGGCCAGTGGATGTGGGTTTCATGAAGGCATTGCATGAGAAGGTCAACATCGTGCCTCT
CATCGCCAAAGCTGACTGTCTTGTCCCCAGTGAGATCCGGAAGCTGAAGGAGCGGATCCGGGAGGAGATT
GACAAGTTTGGGATCCATGTATACCAAGTTCCTGAGTGTGACTCGGACGAGGATGAGGACTTCAAGCAGC
AGGACCGGGAAGTGAAGGAGAGCGCGCCCTTCGCCGTTATAGGCAGCAACACGGTGGTGGAGGCCAAGGG
GCAGCGGGTCCGGGGCCGACTGTACCCCTGGGGATCGTGGAGGTGGAGAACCAGGCGCATTGCGACTTC
GTGAAGCTGCGCAACATGCTCATCCGCACGCATATGCACGACCTCAAGGACGTGACGTGCGACGTGCACT
ACGAGAACTACCGCGCGCACTGCATCCAGCAGATGACCAGCAAACTGACCCAGGACAGCCGCATGGAGAG
CCCCATCCCGATCCTGCCGCTGCCCACCCGACGCGGAGACTGAGAAGCTTATCAGGATGAAGGATGAG
GAACTGAGGCGCATGCAGGAGATGCTGCAGAGGATGAAGCAGCAGATGCAGGACCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206831 protein sequence Red=Cloning site Green=Tags(s) MSTGLRYKSKLATPEDKQDIDKQYVGFATLPNQVHRKSVKKGDFFTLMVAGESGLGKSTLVHSLFLTDLY KDRKLLSAEERISQTV EILKHTVDIEEKGVKLLKLTIVDTPGFGDAVNNTWCWKPITDYVDQQFEQYFRDE SGLNRKNIQDNRVHCCLYFISPFHGHLRPVDVGFMKALHEKVNI VPLIAKADCLVPSEIRKLKERIREEI DKFGIHVYQFPECDSDEDEDFKQDRELKESAPFAVIGSNTVVEAKGQVRGRLYPWGIVEVENQAHCDF VKLRNMLIRTHMDLKDVTCDVHYENYRAHCIQMTSKLTQDSRMESPIPIPLPTDAETEKILRMKDE ELRRMQEMLQRMKQMQDQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002688
ORF Size:	1107 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_002688.4
RefSeq Size:	2090 bp
RefSeq ORF:	1110 bp
Locus ID:	5413
UniProt ID:	Q99719
Cytogenetics:	22q11.21
MW:	42.8 kDa