

Product datasheet for **RC206832A1V**

Human TMEM111 (EMC3) (NM_018447) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TMEM111 (EMC3) (NM_018447) AAV Particle
Tag:	Myc-DDK
Symbol:	TMEM111
Synonyms:	POB; TMEM111
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206832 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCAGGGCCAGAACTGTTGCTCGACTCCAACATCCGCCTCTGGGTGGTCTTACCCATCGTTATCATCA
CTTTCTTCGTAGGCATGATCCGCCACTACGTGTCCATCCTGCTGCAGAGCGACAAGAAGCTACCCAGGA
ACAAGTATCTGACAGTCAAGTCCTAATTCGAAGCAGAGTCCTCAGGGAAAATGGAATAACATTCCTCAAA
CAGTCTTTCTTGACACGAAAATATTATTTCAACAACCCAGAGGATGGATTTTCAAAAACTAAACGGA
AGGTAGTGCCACCTTCTCCTATGACTGATCCTACTATGTTGACAGACATGATGAAAGGGAATGTAACAAA
TGTCCTCCTATGATTCTTATTGGTGGATGGATCAACATGACATTCTCAGGCTTTGTCAACAACCAAGGTC
CCATTTCCACTGACCCTCCGTTTTAAGCCTATGTTACAGCAAGGAATCGAGCTACTCACATTAGATGCAT
CCTGGGTGAGTTCTGCATCCTGGTACTTCTCAATGTATTTGGGCTTCGGAGCATTTACTCTCTGATTCT
GGGCCAAGATAATGCCGCTGACCAATCACGAATGATGCAGGAGCAGATGACGGGAGCAGCCATGGCCATG
CCCGCAGACACAAACAAAGCTTTCAAGACAGAGTGGGAAGCTTTGGAGCTGACGGATCACCAGTGGGCAC
TAGATGATGTGGAAGAAGAGCTCATGGCCAAAGACCTCCACTTCGAAGGCATGTTCAAAAAGGAATTACA
GACCTCTATTTTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206832 protein sequence Red=Cloning site Green=Tags(s) MAGPELLLLDSNIRLWVVLPIVIITFFVGMIRHYVSILLQSDKKLTQEQVSDSQVLIRSRVLRNGKYIPK QSFLTRKYYFNNPEDGFFKKTTRKVVPPSPMTDPTMLTDMMKGNVTNVLPMILIGGWINMTFSGFVTTKV PFPLTLRFKPLMQGIELLTLDASWVSSASWYFLNVFGLRSIYSLILGQDNAADQSRMMQEQMTGAAMAM PADTNKAFKTEWEALELTDHQWALDDVEEELMAKDLHFEGMFKKELQTSIF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_018447
ORF Size:	783 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_018447.1</u>
RefSeq Size:	1107 bp
RefSeq ORF:	786 bp
Locus ID:	55831
UniProt ID:	<u>Q9P0I2</u>
Cytogenetics:	3p25.3
MW:	30 kDa