

Product datasheet for **RC206844L4V**

AIG1 (NM_016108) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AIG1 (NM_016108) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AIG1
Synonyms:	AIG-1; dj95L4.1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016108
ORF Size:	714 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206844).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016108.2
RefSeq Size:	1479 bp
RefSeq ORF:	717 bp
Locus ID:	51390
UniProt ID:	Q9NVV5
Cytogenetics:	6q24.2
Domains:	Far-17a_AIG1
Protein Families:	Transmembrane



[View online »](#)

MW: 27.5 kDa

Gene Summary: May play a role in androgen-regulated growth of hair follicles.[UniProtKB/Swiss-Prot Function]