

Product datasheet for **RC206855L3V**

ACADS (NM_000017) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ACADS (NM_000017) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ACADS
Synonyms:	ACAD3; SCAD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_000017
ORF Size:	1236 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206855).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000017.1
RefSeq Size:	1934 bp
RefSeq ORF:	1239 bp
Locus ID:	35
UniProt ID:	P16219
Cytogenetics:	12q24.31
Domains:	Acyl-CoA_dh, Acyl-CoA_dh_M, Acyl-CoA_dh_N
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Butanoate metabolism, Fatty acid metabolism, Metabolic pathways, Valine, leucine and isoleucine degradation

MW: 44.3 kDa

Gene Summary: This gene encodes a tetrameric mitochondrial flavoprotein, which is a member of the acyl-CoA dehydrogenase family. This enzyme catalyzes the initial step of the mitochondrial fatty acid beta-oxidation pathway. Mutations in this gene have been associated with short-chain acyl-CoA dehydrogenase (SCAD) deficiency. Alternative splicing results in two variants which encode different isoforms. [provided by RefSeq, Oct 2014]