

Product datasheet for **RC206856L1V**

SNRPA1 (NM_003090) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SNRPA1 (NM_003090) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SNRPA1
Synonyms:	Lea1; U2A'
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_003090
ORF Size:	765 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206856).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003090.2
RefSeq Size:	1068 bp
RefSeq ORF:	768 bp
Locus ID:	6627
UniProt ID:	P09661
Cytogenetics:	15q26.3
Domains:	LRR, LRRcap
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Spliceosome

MW: 28.4 kDa

Gene Summary: Involved in pre-mRNA splicing as component of the spliceosome (PubMed:11991638, PubMed:28502770, PubMed:28781166, PubMed:28076346, PubMed:27035939). Associated with sn-RNP U2, where it contributes to the binding of stem loop IV of U2 snRNA (PubMed:9716128, PubMed:27035939).[UniProtKB/Swiss-Prot Function]