

## Product datasheet for **RC206860L2V**

### HPD (NM\_002150) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HPD (NM_002150) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HPD
Synonyms:	4-HPPD; 4HPPD; GLOD3; HPPDASE; PPD
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_002150
ORF Size:	1179 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206860).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002150.2</a> , <a href="#">NP_002141.1</a>
RefSeq Size:	1440 bp
RefSeq ORF:	1182 bp
Locus ID:	3242
UniProt ID:	<a href="#">P32754</a>
Cytogenetics:	12q24.31
Domains:	Glyoxalase
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:** Metabolic pathways, Phenylalanine metabolism, Tyrosine metabolism, Ubiquinone and other terpenoid-quinone biosynthesis

**MW:** 44.9 kDa

**Gene Summary:** The protein encoded by this gene is an enzyme in the catabolic pathway of tyrosine. The encoded protein catalyzes the conversion of 4-hydroxyphenylpyruvate to homogentisate. Defects in this gene are a cause of tyrosinemia type 3 (TYRO3) and hawkinsinuria (HAWK). Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2010]