

Product datasheet for **RC206895A1V**

Human DLX1 (NM_178120) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human DLX1 (NM_178120) AAV Particle
Tag: Myc-DDK
Symbol: DLX1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC206895 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACCATGACCACCATGCCAGAAAGTCTCAACAGCCCCGTGTCGGGCAAGGCGGTGTTTATGGAGTTTG
GGCCGCCCAACCAGCAAATGTCTCTTCTCCATGTCCACGGGCACTACTCCATGCACTGTTTACTCTC
GGCGGGCCATTGCGAGCCCGACGGCGCTACAGCTCAGCCTCGTCTTCTCCGACCGTGGGCTACCCC
TAGTCAACTCGGTGAGCAGCCACGCATCCAGCCCTACATCAGTTCGGTGCAGTCTACCCGGGACGCG
CCAGCCTCGCCAGAGCCGCCTGGAGGACCCAGGGCGGACTCGGAGAAGAGCACGGTGGTGAAGGCGG
TGAAGTGCCTTCAATGGCAAGGAAAAAGATCCGTAACCCAGGACGATTTATTGCAGTTTGCAGTTG
CAGGCTTTGAACCGGAGGTTCCAGCAAACTCAGTACCTAGCTCTGCCGGAGAGGGCGGAGCTCGCGGCCT
CTTTGGGACTCACACAGACTCAGGTCAAGATCTGGTTCCAAAACAAGCGATCCAAGTTCAAGAAGCTGAT
GAAGCAGGGTGGGGCGGCTCTGGAGGGTAGTGCCTGGCCAACGGTCGGGCCCTGTCTGCTGGCTCCCCA
CCCGTGGCCGCCGGCTGGAACCTAACTCTTATCCGGGAAGGGCTCAGGAGGAAACGCGGGCTCCTATA
TCCCCAGCTACACATCGTGGTACCCTTCAGCGCACCAAGAAGCTATGCAGCAACCCCACTTATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC206895 protein sequence
Red=Cloning site Green=Tags(s)

MTMTTPESLNSPVSGKAVFMEFGPPNQMSPPMSHGHYSMHCLHSAGHSQPDGAYSSASSFSRPLGYP
YVNSVSSHASSPYISSVQSYPGSASLAQSRLEDPGADSEKSTVVEGGEVRFNGKGKKIRKPRTIYCSLQL
QALNRRFQQTQYLALPERAELAAASLGLTQTQVKIWFQNKRSKFKKLMKQGGALEGSALANGRALSAGSP
PVPPGWNPNSSSGKSGGNAGSYIPSYTSWYPSAHQEAHQPPQLM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_178120
ORF Size:	765 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_178120.4</u>
RefSeq Size:	2403 bp
RefSeq ORF:	768 bp
Locus ID:	1745
UniProt ID:	<u>P56177</u>
Cytogenetics:	2q31.1
MW:	27.3 kDa