

Product datasheet for **RC206895L2V**

DLX1 (NM_178120) Human Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	DLX1 (NM_178120) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DLX1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_178120
ORF Size:	765 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206895).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178120.4
RefSeq Size:	2403 bp
RefSeq ORF:	768 bp
Locus ID:	1745
UniProt ID:	P56177
Cytogenetics:	2q31.1
Protein Families:	ES Cell Differentiation/IPS, Transcription Factors
MW:	27.3 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of a homeobox transcription factor gene family similar to the *Drosophila* distal-less gene. The encoded protein is localized to the nucleus where it may function as a transcriptional regulator of signals from multiple TGF- β superfamily members. The encoded protein may play a role in the control of craniofacial patterning and the differentiation and survival of inhibitory neurons in the forebrain. This gene is located in a tail-to-tail configuration with another member of the family on the long arm of chromosome 2. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]

Product images:

[RC206895L2] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC206895L2V particle to overexpress human DLX1-mGFP fusion protein.