

Product datasheet for **RC206923L4V**

Mohawk homeobox (MKX) (NM_173576) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Mohawk homeobox (MKX) (NM_173576) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | MKX |
| Synonyms: | C10orf48; IFRX; IRXL1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_173576 |
| ORF Size: | 1056 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC206923). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_173576.1 , NP_775847.1 |
| RefSeq Size: | 3658 bp |
| RefSeq ORF: | 1059 bp |
| Locus ID: | 283078 |
| UniProt ID: | Q8IYA7 |
| Cytogenetics: | 10p12.1 |
| Domains: | homeobox |
| Protein Families: | Transcription Factors |



[View online »](#)

MW: 39.4 kDa

Gene Summary: The protein encoded by this gene is an IRX family-related homeobox protein that may play a role in cell adhesion. Studies in mice suggest that this protein may be a regulator of tendon development. Two transcript variants encoding the same protein have been found for this gene.[provided by RefSeq, Jun 2011]