

Product datasheet for **RC206946A1V**

Human RAB37 (NM_175738) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RAB37 (NM_175738) AAV Particle
Tag: Myc-DDK
Symbol: RAB37
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC206946 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGACCTGCAGAGACCCGATTCTACCAGGGAGGAGCTGGCCCTGACTTCAACGACCACGTCTGCATA
 AGACCATCCTGGTGGGTGACAGTGGTGTGGAAAGACGTCTCTGCTGGTTCAGTTCGATCAGGGCAAGTT
 CATCCCCGGCTCCTTCTCGGCCACTGTGGGCATCGGATTCACGAACAAGGTGGTACTGTGGATGGCGTG
 AGAGTGAAGCTGCAGATCTGGGACACCCTGGGCAGGAACGGTCCGAAGCGTCACCCATGCTTATTACA
 GAGATGCTCAGGCCTTGCTTCTGCTGTATGACATCACCAACAAATCTTCTTCGACAACATCAGGGCCTG
 GCTCACTGAGATTGATGAGTATGCCAGAGGGACGTGGTGATCATGCTGCTAGGCAACAAGCGGATATG
 AGCAGCGAAAGAGTGATCCGTTCCGAAGACGGAGAGACCTTGGCCAGGGAGTACGGTGTTCCTTCCTGG
 AGACCAGCGCCAAGACTGGCATGAATGTGGAGTTAGCCTTTCTGGCCATCGCCAAGGAAGTGAATACCG
 GGCCGGGCATCAGGCGGATGAGCCAGCTTCCAGATCCGAGACTATGTAGAGTCCAGAAGAAGCGCTCC
 AGCTGCTGCTCCTTCATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTAA

Protein Sequence: >RC206946 protein sequence
 Red=Cloning site Green=Tags(s)

MDLQRPDSYQGGAGPDFNDHVLHKTILVGDSGVGKTSLLVQFDQKFIPIGFSFATVGIGFTNKVVTVDGV
 RVKLQIWDTAGQERFRSVTHAYYRDAQALLLYDITNKSSFDNIRAWLTEIHEYAQRDVVIMLLGNKADM
 SSERVIRSEDGETLAREYGVFPFLETSKTMNVELAFLAIAKELKYRAGHQADEPSFQIRDYVESQKKRS
 SCCSFM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



View online »

Serotype:	AAV-2
ACCN:	NM_175738
ORF Size:	648 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_175738.3</u>
RefSeq Size:	3091 bp
RefSeq ORF:	651 bp
Locus ID:	326624
UniProt ID:	<u>Q96AX2</u>
Cytogenetics:	17q25.1
MW:	24.2 kDa