

Product datasheet for **RC206949L3V**

AGBL2 (NM_024783) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | AGBL2 (NM_024783) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | AGBL2 |
| Synonyms: | CCP2 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_024783 |
| ORF Size: | 2706 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC206949). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_024783.2 |
| RefSeq Size: | 3613 bp |
| RefSeq ORF: | 2709 bp |
| Locus ID: | 79841 |
| UniProt ID: | Q5U5Z8 |
| Cytogenetics: | 11p11.2 |
| Protein Families: | Druggable Genome |
| MW: | 104.2 kDa |



[View online »](#)

Gene Summary:

Metallocarboxypeptidase that mediates deglutamylation of target proteins. Catalyzes the deglutamylation of polyglutamate side chains generated by post-translational polyglutamylation in proteins such as tubulins. Also removes gene-encoded polyglutamates from the carboxy-terminus of target proteins such as MYLK. Does not show detyrosinase or deglycylase activities from the carboxy-terminus of tubulin.[UniProtKB/Swiss-Prot Function]