

Product datasheet for **RC206962A1V**

Human IL17D (NM_138284) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human IL17D (NM_138284) AAV Particle
Tag: Myc-DDK
Symbol: IL17D
Synonyms: IL-17D
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC206962 representing NM_138284
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCTGGTAGCCGGCTTCCTGCTGGCGCTGCCGCCGAGCTGGGCCGCGGCGCCCCGAGGGCGGGCAGGC
GCCCCGCGCGGCCGCGGGGCTGCGCGGACCGCCGAGGAGCTACTGGAGCAGCTGTACGGGCGCCTGGC
GGCCGGCGTGCTCAGTGCCTTCCACCACACGCTGCAGCTGGGGCCGCGTGAGCAGGCGCGCAACGCGAGC
TGCCCGGCAGGGGGCAGGCCCCGCCGACCGCCGCTTCCGGCCGCCACCAACCTGCGCAGCGTGTGCGCCT
GGGCCTACAGAATCTCCTACGACCGGGCGAGGTACCCAGGTACCTGCCTGAAGCCTACTGCCTGTGCCG
GGGCTGCCTGACCGGGCTGTTCCGCGAGGAGGACGTGCGCTTCCGACGCGCCCTGTCTACATGCCCAACC
GTGCTCCTGCGCCGACCCCGCCTGCGCCGGCGGCGTTCGCTCTACCCGAGGCCTACGTACCATCC
CCGTGGGCTGCACCTGCGTCCCCGAGCCGAGAGGACGACAGCAGCATCAACTCCAGCATCGACAAACA
GGGCGCCAAGCTCCTGCTGGGCCCAACGACGCGCCCGCTGGCCCC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC206962 representing NM_138284
Red=Cloning site Green=Tags(s)

MLVAGFLLALPPSWAAGAPRAGRRPARPRGCADRPEELLEQLYGRLAAGVLSAFHHTLQLGPREQARNAS
CPAGGRPADDRFRPPTNLRVSPWAYRISYDPARYPRYLPEAYCLCRGCLTGLFGEDVRFRRSAPVYMP
VVLRRTPACAGGRSVYTEAYVTIPVGCTCVPEPEKDADSINSSIDKQAKLLLPNDAPAGP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_138284
ORF Size:	606 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_138284.1</u> , <u>NP_612141.1</u>
RefSeq Size:	1873 bp
RefSeq ORF:	609 bp
Locus ID:	53342
UniProt ID:	<u>Q8TAD2</u>
Cytogenetics:	13q12.11
MW:	21.4 kDa