

Product datasheet for **RC206991A1V**

Human CD94 (KLRD1) (NM_002262) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CD94 (KLRD1) (NM_002262) AAV Particle
Tag: Myc-DDK
Symbol: CD94
Synonyms: CD94
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC206991 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCAGTGTTTAAGACCACTCTGTGGAGGTTAATTTCTGGGACCTTAGGGATAATATGCCTTTCGTTGA
 TGGCTACGTTGGGAATTTGTTGAAAAATCTTTTACTAACTGAGTATTGAGCCAGCATTACTCCAGG
 ACCCAACATAGAACTCCAGAAAGACTCTGACTGCTGTTCTTGCCAAGAAAAATGGGTTGGGTACCGGTGC
 AACTGTTACTTCATTTCCAGTGAACAGAAAACTTGAACGAAAGTCGGCATCTCTGTGCTTCTCAGAAAT
 CCAGCCTGCTTCAGCTTCAAAACACAGATGAAGTGGATTTATGAGCTCCAGTCAACAATTTTACTGGAT
 TGGACTCTTTACAGTGAGGAGCACACCGCTGGTTGTGGGAGAATGGCTCTGCACTCTCCAGTATCTA
 TTTCCATCATTTGAACTTTTAATACAAAGAACTGCATAGCGTATAATCCAATGGAAATGCTTTAGATG
 AATCCTGTGAAGATAAAATCGTTATATCTGTAAGCAACAGCTCATT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC206991 protein sequence
 Red=Cloning site Green=Tags(s)

MAVFKTTLWRLISGTLGIICLSMATLGILLKNSFTKLSIEPAFTPGPNIELQKDSGCCSCQEKWVGRC
 NCYFISSEQKTWNEsrhlcasqkssllqlqntdelDFMSSSQFYWIGLSYSEEHTAWLWENGSAISQYL
 FPSFETFTKNCIAYNPNGNALDESCEDKNRYICKQQLI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_002262
ORF Size:	537 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_002262.2</u>
RefSeq Size:	3258 bp
RefSeq ORF:	540 bp
Locus ID:	3824
UniProt ID:	<u>Q13241</u>
Cytogenetics:	12p13.2
MW:	20.5 kDa