

Product datasheet for **RC207003L4V**

LILRA3 (NM_006865) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LILRA3 (NM_006865) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LILRA3
Synonyms:	CD85E; HM31; HM43; ILT-6; ILT6; LIR-4; LIR4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006865
ORF Size:	1317 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207003).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006865.2 , NP_006856.2
RefSeq Size:	1608 bp
RefSeq ORF:	1320 bp
Locus ID:	11026
UniProt ID:	Q8N6C8
Cytogenetics:	19q13.4
Domains:	ig, IGc2, IG
Protein Families:	Secreted Protein



[View online »](#)

MW: 47.5 kDa

Gene Summary: This gene encodes a member of a family of immunoreceptors that are expressed predominantly in monocytes and B cells, and at lower levels in dendritic cells and natural killer cells. The encoded protein lacks the transmembrane region found in other members of this family. It acts as a soluble receptor for class I major histocompatibility complex (MHC) antigens. Alternatively spliced transcript variants encoding different isoforms have been found. This gene is located in a cluster of related genes on chromosome 19 and is polymorphic in human populations, with many individuals containing a deletion of this genomic region. [provided by RefSeq, Mar 2014]