

## Product datasheet for **RC207004L3V**

### Fos B (FOSB) (NM\_006732) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Fos B (FOSB) (NM_006732) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Fos B
Synonyms:	AP-1; G0S3; GOS3; GOSB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006732
ORF Size:	1014 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207004).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006732.1</a>
RefSeq Size:	3776 bp
RefSeq ORF:	1017 bp
Locus ID:	2354
UniProt ID:	<a href="#">P53539</a>
Cytogenetics:	19q13.32
Domains:	BRLZ
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**MW:** 35.9 kDa

**Gene Summary:** The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]