

### **Product datasheet for RC207015L2**

# NR2E1 (NM\_003269) Human Tagged Lenti ORF Clone

#### **Product data:**

**Product Type:** Expression Plasmids

**Product Name:** NR2E1 (NM\_003269) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: NR2E1

Synonyms: TLL; TLX; XTLL

Mammalian Cell None

Selection:

**Vector:** pLenti-C-mGFP (PS100071)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC207015).

Sequence:

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 





st The last codon before the Stop codon of the ORF.

**ACCN:** NM\_003269

ORF Size: 1155 bp



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



#### NR2E1 (NM\_003269) Human Tagged Lenti ORF Clone - RC207015L2

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 003269.2</u>

 RefSeq Size:
 3225 bp

 RefSeq ORF:
 1158 bp

 Locus ID:
 7101

UniProt ID: Q9Y466

Cytogenetics: 6q21

**Domains:** HOLl, zf-C4

**Protein Families:** Druggable Genome, Nuclear Hormone Receptor, Transcription Factors

MW: 42.4 kDa


**Gene Summary:** The protein encoded by this gene is an orphan receptor involved in retinal development. The

encoded protein also regulates adult neural stem cell proliferation and may be involved in control of aggressive behavior. Two transcript variants encoding different isoforms have been

found for this gene. [provided by RefSeq, Aug 2015]



## **Product images:**



Circular map for RC207015L2