

Product datasheet for **RC207018A1V**

Human CXCL7 (PPBP) (NM_002704) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CXCL7 (PPBP) (NM_002704) AAV Particle
Tag:	Myc-DDK
Symbol:	CXCL7
Synonyms:	B-TG1; Beta-TG; CTAP-III; CTAP3; CTAPIII; CXCL7; LA-PF4; LDGF; MDGF; NAP-2; PBP; SCYB7; TC1; TC2; TGB; TGB1; THBGB; THBGB1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207018 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC ATGAGCCTCAGACTTGATACCACCCCTTCCTGTAAAGTGCAGAGCACTTCATGCCTTGACAGGTGCTGCTGCTTCTGCTCATTGCTGCTGACTGCTCTGGCTTCCTCCACCAAAGGACAACTAAGAGAACTTGCGCAAAGGCAAAGAGGAAAGTCTAGACAGTGACTTGATGCTGAACTCCGCTGCATGTGTATAAAGACAACCTCTGGAATTCATCCCAAAACATCCAAAGTTTGAAGTGATCGGGAAGGAACCCATTGCAACCAAGTCGAAGTGATAGCCACACTGAAGGATGGGAGGAAATCTGCCTGGACCCAGATGCTCCAGAATCAAGAAAATTGTACAGAAAAATTGGCAGGTGATGAATCTGCTGAT ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC207018 protein sequence Red=Cloning site Green=Tags(s) MSLRDLTPSCNSARPLHALQVLLLLSLLLTALASSTKGQTKRNLAKGKEESLSDLYAELRCMCIKTTSGIHPKNIQSLEVIGKTHCNQVEVIATLKDGRKICLDPDAPRIKKIVQKKLAGDESAD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002704


[View online »](#)

ORF Size:	384 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_002704.2</u>
RefSeq Size:	1307 bp
RefSeq ORF:	387 bp
Locus ID:	5473
UniProt ID:	<u>P02775</u>
Cytogenetics:	4q13.3
MW:	13.9 kDa