

Product datasheet for **RC207021L3V**

Hexokinase Type III (HK3) (NM_002115) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hexokinase Type III (HK3) (NM_002115) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Hexokinase Type III
Synonyms:	HKIII; HXK3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002115
ORF Size:	2769 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207021).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002115.1
RefSeq Size:	3062 bp
RefSeq ORF:	2772 bp
Locus ID:	3101
UniProt ID:	P52790
Cytogenetics:	5q35.2
Domains:	hexokinase
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways:	Amino sugar and nucleotide sugar metabolism, Fructose and mannose metabolism, Galactose metabolism, Glycolysis / Gluconeogenesis, Insulin signaling pathway, Metabolic pathways, Starch and sucrose metabolism, Type II diabetes mellitus
MW:	98.8 kDa
Gene Summary:	Hexokinases phosphorylate glucose to produce glucose-6-phosphate, the first step in most glucose metabolism pathways. This gene encodes hexokinase 3. Similar to hexokinases 1 and 2, this allosteric enzyme is inhibited by its product glucose-6-phosphate. [provided by RefSeq, Apr 2009]