

Product datasheet for RC207023L2

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

CD33 (NM_001772) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: CD33 (NM_001772) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: CD33

Synonyms: p67; SIGLEC-3; SIGLEC3

Mammalian Cell None

Selection:

Vector: pLenti-C-mGFP (PS100071)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC207023).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_001772

ORF Size: 1092 bp



CD33 (NM_001772) Human Tagged Lenti ORF Clone - RC207023L2

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

> reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: NM 001772.2, NP 001763.2

945

RefSeq Size: 1466 bp RefSeq ORF: 1095 bp

UniProt ID: P20138

Locus ID:

Cytogenetics:

19q13.41 Druggable Genome, Transmembrane **Protein Families:**

Protein Pathways: Hematopoietic cell lineage

MW: 39.7 kDa

Gene Summary: Sialic-acid-binding immunoglobulin-like lectin (Siglec) that plays a role in mediating cell-cell

interactions and in maintaining immune cells in a resting state (PubMed:10611343,

PubMed:15597323, PubMed:11320212). Preferentially recognizes and binds alpha-2,3- and more avidly alpha-2,6-linked sialic acid-bearing glycans (PubMed:7718872). Upon engagement

of ligands such as C1q or syalylated glycoproteins, two immunoreceptor tyrosine-based inhibitory motifs (ITIMs) located in CD33 cytoplasmic tail are phosphorylated by Src-like

kinases such as LCK (PubMed:28325905, PubMed:10887109). These phosphorylations provide


docking sites for the recruitment and activation of protein-tyrosine phosphatases

PTPN6/SHP-1 and PTPN11/SHP-2 (PubMed:10556798, PubMed:10206955, PubMed:10887109). In turn, these phosphatases regulate downstream pathways through dephosphorylation of signaling molecules (PubMed:10206955, PubMed:10887109). One of the repressive effect of CD33 on monocyte activation requires phosphoinositide 3-kinase/PI3K (PubMed:15597323).


[UniProtKB/Swiss-Prot Function]



Product images:



Circular map for RC207023L2



Double digestion of RC207023L2 using Sgfl and Mlul $\,$