

Product datasheet for **RC207037L3V**

Alpha SNAP (NAPA) (NM_003827) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Alpha SNAP (NAPA) (NM_003827) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | Alpha SNAP |
| Synonyms: | SNAPA |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_003827 |
| ORF Size: | 885 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC207037). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_003827.2 |
| RefSeq Size: | 1865 bp |
| RefSeq ORF: | 888 bp |
| Locus ID: | 8775 |
| UniProt ID: | P54920 |
| Cytogenetics: | 19q13.32-q13.33 |
| Domains: | NSF |
| MW: | 33.2 kDa |


[View online »](#)

Gene Summary:

This gene encodes a member of the soluble NSF attachment protein (SNAP) family. SNAP proteins play a critical role in the docking and fusion of vesicles to target membranes as part of the 20S NSF-SNAP-SNARE complex. The encoded protein plays a role in the completion of membrane fusion by mediating the interaction of N-ethylmaleimide-sensitive factor (NSF) with the vesicle-associated and membrane-associated SNAP receptor (SNARE) complex, and stimulating the ATPase activity of NSF. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jun 2011]