

Product datasheet for **RC207051L4V**

Kir2.2 (KCNJ12) (NM_021012) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kir2.2 (KCNJ12) (NM_021012) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Kir2.2
Synonyms:	hIRK; hIRK1; hkir2.2x; IRK-2; IRK2; kcnj12x; KCNJN1; Kir2.2; Kir2.2v
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021012
ORF Size:	1299 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207051).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021012.4
RefSeq Size:	5230 bp
RefSeq ORF:	1302 bp
Locus ID:	3768
UniProt ID:	Q14500
Cytogenetics:	17p11.2
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane
MW:	48.8 kDa



[View online »](#)

Gene Summary:

This gene encodes an inwardly rectifying K⁺ channel which may be blocked by divalent cations. This protein is thought to be one of multiple inwardly rectifying channels which contribute to the cardiac inward rectifier current (IK1). The gene is located within the Smith-Magenis syndrome region on chromosome 17. [provided by RefSeq, Jul 2008]