

## Product datasheet for **RC207061L3V**

### SLC25A34 (NM\_207348) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC25A34 (NM_207348) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC25A34
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_207348
ORF Size:	912 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207061).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_207348.1</a>
RefSeq Size:	3172 bp
RefSeq ORF:	915 bp
Locus ID:	284723
UniProt ID:	<a href="#">Q6PIV7</a>
Cytogenetics:	1p36.21
Protein Families:	Transmembrane
MW:	32.2 kDa



[View online »](#)

**Gene Summary:** SLC25A34 belongs to the SLC25 family of mitochondrial carrier proteins (Haitina et al., 2006 [PubMed 16949250]).[supplied by OMIM, Mar 2008]