

Product datasheet for **RC207067L3V**

TMEM139 (NM_153345) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TMEM139 (NM_153345) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM139
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153345
ORF Size:	648 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207067).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153345.1
RefSeq Size:	1980 bp
RefSeq ORF:	651 bp
Locus ID:	135932
UniProt ID:	Q8IV31
Cytogenetics:	7q34
Protein Families:	Transmembrane
MW:	23.7 kDa



[View online »](#)

Gene Summary:

May be involved in cellular trafficking of proteins such as SLC4A1.[UniProtKB/Swiss-Prot Function]