

## Product datasheet for **RC207075L3V**

### MRPS11 (NM\_022839) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | MRPS11 (NM_022839) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | MRPS11   |
| Synonyms:                 | HCC-2; MRP-S11; S11mt  |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_022839  |
| ORF Size:                 | 582 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC207075).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_022839.2</a>  |
| RefSeq Size:              | 1501 bp  |
| RefSeq ORF:               | 585 bp   |
| Locus ID:                 | 64963  |
| UniProt ID:               | <a href="#">P82912</a>   |
| Cytogenetics:             | 15q25.3  |
| Domains:                  | Ribosomal_S11  |
| MW:                       | 20.6 kDa   |



[View online »](#)

**Gene Summary:**

Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein that contains a high level of sequence similarity with ribosomal protein S11P family members. A pseudogene corresponding to this gene is found on chromosome 20. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2016]