

Product datasheet for **RC207080L1V**

Myoneurin (MYNN) (NM_018657) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Myoneurin (MYNN) (NM_018657) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Myoneurin
Synonyms:	OSZF; SBBIZ1; ZBTB31; ZNF902
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_018657
ORF Size:	1830 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207080).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018657.3
RefSeq Size:	5300 bp
RefSeq ORF:	1833 bp
Locus ID:	55892
UniProt ID:	Q9NPC7
Cytogenetics:	3q26.2
Domains:	BTB, zf-C2H2
Protein Families:	Transcription Factors



[View online »](#)

MW: 68.7 kDa

Gene Summary: This gene encodes a member of the BTB/POZ and zinc finger domain-containing protein family that are involved in the control of gene expression. Alternative splicing results in multiple transcript variants and a pseudogene has been identified on chromosome 14. [provided by RefSeq, Jun 2010]