

Product datasheet for **RC207089A1V**

Human OSGEP (NM_017807) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human OSGEP (NM_017807) AAV Particle
Tag:	Myc-DDK
Symbol:	OSGEP
Synonyms:	GAMOS3; GCPL1; KAE1; OSGEP1; PRSMG1; TCS3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207089 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCCGGCGGTGCTGGGTTTTGAAGGCAGCGCCAATAAGATTGGCGTGGGCGTGGTGCGGGATGGCAAGG
TGCTGGCGAACC CGCGCGGACTTACGTACGCCTCCTGGCACAGGATTCCTTCCAGGTGATACAGCCAG
GCATACCCGAGCTGTTATCCTAGACCTGCTGCAGGAGGCACTAACAGAGTCTGGATTAACCTCCCAGGAT
ATCGACTGCATTGCATACCAAGGGCCCTGGCATGGGTGCCCACTGGTTTCTGTGGCTGTTGTGGCCC
GTACTGTGGCCCACTGTGGAATAAGCCATTGGTGGGTGTGAACCACTGTATAGGCCACATTGAGATGGG
CCGCCTCATCACTGGAGCCACCAGCCCAACCGTGTGTATGTGAGTGGAGGAAATACGCAGGTGATTGCA
TACTCGGAACATCGTTACCGTATCTTTGGGGAACCATCGATATTGCAGTGGTAATTGTCTGGATCGTT
TTGCTCGAGTGTGAAGATTTCTAACGACCAAGTCCAGGATACAACATTGAACAGATGGCAAAGCGAGG
CAAGAAGCTAGTTGAGCTGCCATACCTGTAAAGGGGATGGACGTCTATTCTCAGGGATCCTGTCTTTC
ATTGAGGATGTAGCCATCGGATGTGGCCACAGGCGAGTGTACTCTGAGGATCTGTGTTTCTCCCTGC
AGGAACTGTGTTTGAATGCTGGTAGAGATCACAGAGCGAGCCATGGCACATTGTGGCTCCCAGGAGGC
CCTCATTGTGGGAGGAGTGGGTGTAATGTGAGGCTACAGGAGATGATGGCAACAATGTGCCAGGAACGT
GGAGCCCGGCTTTTGTACAGATGAGAGATTCTGTATTGACAATGGAGCGATGATAGCCAGGCTGGCT
GGGAGATGTTTTCGGGCTGGACACAGGACCCCACTCAGTGATTCTGGGGTTACACAGAGGTATCGGACAGA
TGAAGTAGAGGTGACCTGGAGGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207089 protein sequence Red=Cloning site Green=Tags(s) MPAVLGFEGSANKIGVGVVRDQKVLNPRRTYVTPPGTGFLPGDTARHHRVILDLLQEALTESGLTSQD IDCIAYTKGPGMGAPLVSVAVVARTVAQLWNKPLVGVNHCIGHIEMGRLITGATSPTVLYVSGGNTQVIA YSEHRYRIFGETIDIAVGNCLDRFARVLKISNDPSPGYNIEQMAKRGKKLVLPYTVKGMDVSFSGILSF IEDVAHRMLATGECTPEDLCFSLQETVFAMLVEITERAMAHCGSQEALIVGGVGCNVRQLQEMMATMCQER GARLFATDERFCIDNGAMIAQAGWEMFRAGHRTPLSDSGVTQRYRTDEVEVTWRD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_017807
ORF Size:	1005 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_017807.2</u>
RefSeq Size:	1642 bp
RefSeq ORF:	1008 bp
Locus ID:	55644
UniProt ID:	<u>Q9NPF4</u>
Cytogenetics:	14q11.2
MW:	36.4 kDa