

Product datasheet for RC207100L3

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Dishevelled 3 (DVL3) (NM_004423) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Dishevelled 3 (DVL3) (NM_004423) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: Dishevelled 3

Synonyms: DRS3

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC207100).

Sequence:

Sgfl-Mlul

Restriction Sites: Cloning Scheme:






^{*} The last codon before the Stop codon of the ORF.





Plasmid Map:



ACCN: NM_004423
ORF Size: 2148 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 004423.3</u>

RefSeq Size: 5062 bp RefSeq ORF: 2151 bp



Dishevelled 3 (DVL3) (NM_004423) Human Tagged Lenti ORF Clone - RC207100L3

 Locus ID:
 1857

 UniProt ID:
 Q92997

 Cytogenetics:
 3q27.1

Domains: DEP, DAX, PDZ, Dishevelled

Protein Families: Druggable Genome, ES Cell Differentiation/IPS

Protein Pathways: Basal cell carcinoma, Colorectal cancer, Melanogenesis, Notch signaling pathway, Pathways in

cancer, Wnt signaling pathway

MW: 77.9 kDa

Gene Summary: This gene is a member of a multi-gene family which shares strong similarity with the

Drosophila dishevelled gene, dsh. The Drosophila dishevelled gene encodes a cytoplasmic

phosphoprotein that regulates cell proliferation. [provided by RefSeq, Jul 2008]