

Product datasheet for **RC207109A1V**

Human CSRP1 (NM_004078) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CSRP1 (NM_004078) AAV Particle
Tag: Myc-DDK
Symbol: CSRP1
Synonyms: CRP; CRP1; CSRP; CYRP; D1S181E; HEL-141; HEL-S-286
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC207109 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCCGAAGTGGGGAGGAGGCAAGAAATGTGGGGTGTGTCAGAAGACGGTTTACTTTGCCGAAGAGGTTTC
AGTGCGAAGGCAACAGCTTCCATAAATCCTGCTTCCTGTGCATGGTCTGCAAGAAGAATCTGGACAGTAC
CACTGTGGCCGTGCATGGTGAGGAGATTACTGCAAGTCCTGCTACGGCAAGAAGTATGGGCCCAAAGGC
TATGGCTACGGGCAGGGCGCAGGCACCCCTCAGCACTGACAAGGGGGAGTCGCTGGGTATCAAGCACGAGG
AAGCCCCTGGCCACAGGCCACCACCAACCCCAATGCATCCAAATTTGCCAGAAGATTGGTGGCTCCGA
GCGCTGCCCCCGATGCAGCCAGGCAGTCTATGCTGCGGAGAAGGTGATTGGTGCTGGGAAGTCCTGGCAT
AAGGCCTGCTTTTCGATGTGCCAAGTGTGGCAAAGGCCTTGAGTCAACCACCCTGGCAGACAAGGATGGCG
AGATTTACTGCAAAGGATGTTATGCTAAAACTTCGGGCCCAAGGGCTTTGGTTTGGGCAAGGAGCTGG
GGCCTTGGTCCACTCTGAG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC207109 protein sequence
Red=Cloning site Green=Tags(s)

MPNWGGGKKCGVCQKTVYFAEEVQCEGNSFHKSCFLCMVCKKNLDSTTVAVHGEEIYCKSCYKKGYPKG
YGYGQAGTLSTDKGESLGIKHEEAPGHRPTTNPNAKFAQKIGGSERCPSCQAVYAAEKVIGAGKSWH
KACFRCAKCGKLESTTLADKDGEIYCKGKYAKNFGPKGFGFGGAGALVHSE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



Serotype:	AAV-2
ACCN:	NM_004078
ORF Size:	579 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_004078.1</u>
RefSeq Size:	1970 bp
RefSeq ORF:	582 bp
Locus ID:	1465
UniProt ID:	<u>P21291</u>
Cytogenetics:	1q32.1
MW:	20.6 kDa