

Product datasheet for **RC207110L4V**

ACYP1 (NM_001107) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ACYP1 (NM_001107) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ACYP1
Synonyms:	ACYPE
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001107
ORF Size:	297 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207110).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001107.3
RefSeq Size:	621 bp
RefSeq ORF:	300 bp
Locus ID:	97
UniProt ID:	P07311
Cytogenetics:	14q24.3
Domains:	Acylphosphatase
Protein Pathways:	Pyruvate metabolism



[View online »](#)

MW: 11.3 kDa

Gene Summary: This gene is a member of the acylphosphatase family. The encoded protein is a small cytosolic enzyme that catalyzes the hydrolysis of the carboxyl-phosphate bond of acylphosphates. Two isoenzymes have been isolated and described based on their tissue localization: erythrocyte (common) type acylphosphatase encoded by this gene, and muscle type acylphosphatase. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2014]