

## Product datasheet for **RC207111L4V**

### AADAC (NM\_001086) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	AADAC (NM_001086) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AADAC
Synonyms:	CES5A1; DAC
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001086
ORF Size:	1197 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207111).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001086.2</a>
RefSeq Size:	1725 bp
RefSeq ORF:	1200 bp
Locus ID:	13
UniProt ID:	<a href="#">P22760</a>
Cytogenetics:	3q25.1
Protein Families:	Druggable Genome, Transmembrane
MW:	45.7 kDa



[View online »](#)

**Gene Summary:**

Microsomal arylacetamide deacetylase competes against the activity of cytosolic arylamine N-acetyltransferase, which catalyzes one of the initial biotransformation pathways for arylamine and heterocyclic amine carcinogens [provided by RefSeq, Jul 2008]