

Product datasheet for RC207120L4

FHIT (NM_002012) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: FHIT (NM 002012) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: FHIT

Synonyms: AP3Aase; FRA3B

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC207120).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_002012

ORF Size: 441 bp



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



FHIT (NM_002012) Human Tagged Lenti ORF Clone - RC207120L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 002012.1</u>

RefSeq Size: 1103 bp RefSeq ORF: 444 bp Locus ID: 2272

UniProt ID: P49789

Cytogenetics: 3p14.2

Domains: HIT

Protein Pathways: Non-small cell lung cancer, Purine metabolism, Small cell lung cancer

MW: 16.9 kDa

Gene Summary: The protein encoded by this gene is a P1-P3-bis(5'-adenosyl) triphosphate hydrolase involved

in purine metabolism. This gene encompasses the common fragile site FRA3B on

chromosome 3, where carcinogen-induced damage can lead to translocations and aberrant transcripts. In fact, aberrant transcripts from this gene have been found in about half of all


esophageal, stomach, and colon carcinomas. The encoded protein is also a tumor

suppressor, as loss of its activity results in replication stress and DNA damage. [provided by

RefSeq, Aug 2017]



Product images:



Circular map for RC207120L4