

Product datasheet for **RC207126A1V**

Human LYNX1-SLURP2 (NM_023946) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human LYNX1-SLURP2 (NM_023946) AAV Particle
Tag: Myc-DDK
Symbol: LYNX1-SLURP2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC207126 representing NM_023946
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACGCCCTGCTCACCTGATCCTGGTGGTCTCATGGGCTTACCTCTGGCCAGGCCTTGGACTGCC
ACGTGTGTGCCTACAACGGAGACAAGTCTCAACCCCATGCGCTGCCCGGCTATGGTTGCCTACTGCAT
GACCACGCGCACCTCTGCAGCCGAAGCCATATGGTGTACCAGTGCACGGGCTTCGGAGGGTCTCCCAT
GGATCCAGATGCCTGAGGGACTCCACCCACTGTGTACCAGTGCACCCGGGCTCCTCAGAACACCGAGG
ATTTGCCTCTGGTCACCAAGATGTGCCACATAGGCTGCCCGGATATCCCAGCCTGGGCCTGGGCCCTA
CGTATCCATCGTTGCTGCCAGACCAGCCTCTGCAACCATGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC207126 representing NM_023946
Red=Cloning site Green=Tags(s)

MTPLLTLILVVLMLPLAQAALDCHVCAYNQDNCNPNMRCMPAMVAYCMTTRTSAEAIWCHQCTGFGGCSH
GSRCLRSTHCVTATRVLSNTEDLPLVTKMCHIGCPDIPSLGLGPYVSIACQTSKLNHD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_023946
ORF Size: 393 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_023946.1</u>
RefSeq Size:	1352 bp
RefSeq ORF:	396 bp
Locus ID:	111188157
UniProt ID:	<u>PODP58</u>
Cytogenetics:	8q24.3
MW:	13.8 kDa