

Product datasheet for **RC207126L4V**

LYNX1 (LYNX1-SLURP2) (NM_023946) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LYNX1 (LYNX1-SLURP2) (NM_023946) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LYNX1-SLURP2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_023946
ORF Size:	393 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207126).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023946.1
RefSeq Size:	1352 bp
RefSeq ORF:	396 bp
Locus ID:	111188157
UniProt ID:	P0DP58
Cytogenetics:	8q24.3
MW:	13.8 kDa



[View online »](#)

Gene Summary:

This locus represents naturally occurring read-through transcription between the neighboring LYNX1 and SLURP2 genes. The readthrough transcript encodes a fusion protein comprised of sequence sharing identity with each individual gene product. The significance of this read-through transcription and the function of the resulting protein product have not yet been determined. [provided by RefSeq, Sep 2017]