

## Product datasheet for **RC207135L4V**

### PHF17 (JADE1) (NM\_024900) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PHF17 (JADE1) (NM_024900) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PHF17
Synonyms:	PHF17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024900
ORF Size:	1527 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207135).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024900.3</a>
RefSeq Size:	4204 bp
RefSeq ORF:	1530 bp
Locus ID:	79960
UniProt ID:	<a href="#">Q6IE81</a>
Cytogenetics:	4q28.2
Domains:	PHD
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**MW:** 58.4 kDa

**Gene Summary:** Component of the HBO1 complex which has a histone H4-specific acetyltransferase activity, a reduced activity toward histone H3 and is responsible for the bulk of histone H4 acetylation in vivo. Transcriptional coactivator, it may also promote acetylation of nucleosomal histone H4 by KAT5. Promotes apoptosis. May act as a renal tumor suppressor. Negatively regulates canonical Wnt signaling; at least in part, cooperates with NPHP4 in this function.  
[UniProtKB/Swiss-Prot Function]