

Product datasheet for **RC207152L3V**

Cryptochrome I (CRY1) (NM_004075) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cryptochrome I (CRY1) (NM_004075) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Cryptochrome I
Synonyms:	DSPD; PHLL1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004075
ORF Size:	1758 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207152).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004075.2
RefSeq Size:	3310 bp
RefSeq ORF:	1761 bp
Locus ID:	1407
UniProt ID:	Q16526
Cytogenetics:	12q23.3
Domains:	FAD_binding_7, DNA_photolyase
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Circadian rhythm - mammal

MW: 66.4 kDa

Gene Summary: This gene encodes a flavin adenine dinucleotide-binding protein that is a key component of the circadian core oscillator complex, which regulates the circadian clock. This gene is upregulated by CLOCK/ARNTL heterodimers but then represses this upregulation in a feedback loop using PER/CRY heterodimers to interact with CLOCK/ARNTL. Polymorphisms in this gene have been associated with altered sleep patterns. The encoded protein is widely conserved across plants and animals. Loss of the related gene in mouse results in a shortened circadian cycle in complete darkness. [provided by RefSeq, Jan 2014]