

Product datasheet for **RC207162A1V**

Human TIRAP (NM_001039661) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TIRAP (NM_001039661) AAV Particle
Tag:	Myc-DDK
Symbol:	TIRAP
Synonyms:	BACTS1; Mal; MyD88-2; wyatt
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207162 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCATCATCGACCTCCCTCCAGCTCCTGGCTCTCGGCCTAAGAAGCCTCTAGGCAAGATGGCTGACT
GGTTCAGGCAGACCTGCTGAAGAAGCCCAAGAAGAGGCCCAACTCCCCAGAAAGCACCTCCAGCGATGC
TTCACAGCCTACCTCACAGGACAGCCCACTACCCCAAGCCTCAGCTCAGTCACGTCTCCAGCCTGCCA
CCCACACATGCGAGTGACAGTGGCAGTAGTCGCTGGAGCAAAGACTATGACGTCTGCGTGTGCCACAGTG
AGGAAGACCTGGTGGCCGCCAGGACCTGGTCTCCTACTTGGAAGGCAGCACTGCCAGCCTGCGCTGCTT
CCTGCAACTCCGGGATGCAACCCAGGCGGCGCTATAGTGTCCGAGCTGTGCCAGGCACTGAGCAGTAGT
CACTGCCGGGTGCTGCTCATCACGCCGGGCTTCCTTCAGGACCCCTGGTGCAAGTACCAGATGCTGCAGG
CCCTGACCGAGGCTCCAGGGGCCGAGGGCTGCACCATCCCCCTGCTGTCGGGCCTCAGCAGAGCTGCCTA
CCCACCTGAGCTCCGATTATGTACTACGTGATGGCAGGGGCCCTGATGGTGGCTTTCGTCAAGTCAAA
GAAGCTGTATGCGTTATCTGCAGACACTCAGTTGGCACTTGTTATATCATGGGACCCCGGAAATTGGAG
TGAAGCTAGAAACAGAAAACCCATGCAGGGCCTCGGATCCCACAAATGTGACAAGAGGTATAGGGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207162 protein sequence Red=Cloning site Green=Tags(s) MASSTSLPAPGSRPKKPLGKMADWFRQTLLKKPKKRPNSPESTSSDASQPTSQDSPLPPSLSSVTSPSLP PTHASDSGSSRWSKDYDVCVCHSEEDLVAAQDLVSYLEGSTASLRCLQLRDATPGGAI VSEL CQALSSS HCRVLLITPGFLQDPWCKYQMLQALTEAPGAEGCTIPLL SGLSRAAYPPELRFMYVVDGRGPDGGFRQVK EAVMRYLQTL SWHLLYHGTPEIGVKLETENPCRASDSHKCDKRYRE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001039661
ORF Size:	768 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001039661.1</u>
RefSeq Size:	2348 bp
RefSeq ORF:	666 bp
Locus ID:	114609
UniProt ID:	<u>P58753</u>
Cytogenetics:	11q24.2
MW:	28 kDa