

## Product datasheet for **RC207166L4V**

### SPIC (NM\_152323) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SPIC (NM_152323) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SPIC
Synonyms:	SPI-C
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_152323
ORF Size:	744 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207166).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_152323.1</a>
RefSeq Size:	1116 bp
RefSeq ORF:	747 bp
Locus ID:	121599
UniProt ID:	<a href="#">Q8N5J4</a>
Cytogenetics:	12q23.2
Protein Families:	Transcription Factors
MW:	29.2 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene regulates the development of red pulp macrophages, which are necessary for iron homeostasis and the recycling of red blood cells. [provided by RefSeq, Aug 2016]