

Product datasheet for **RC207169A1V**

Human EAF1 (NM_033083) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human EAF1 (NM_033083) AAV Particle
Tag: Myc-DDK
Symbol: EAF1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC207169 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAATGGGACCGCAAACCCGCTGCTGGACCGCGAGGAACATTGCCTGAGGCTCGGGGAGAGCTTCGAGA
AGCGGCCCGGGCCTCTTCCACACTATTCGTTATGATTTTAAACCAGCATCTATAGACACTTCCTGTGA
AGGAGAGCTTCAAGTTGGCAAAGGAGATGAAGTCACAATTACACTTCACATATCCCTGGATCCACACCA
CCCATGACTGTGTTCAAGGGGAACAAACGGCCTTACCAGAAAGACTGTGTGCTTATTATTAATCATGACA
CTGGTGAATTTGTGCTGGAAAACTCAGTAGCAGCATTAGGTGAAGAAAACAAGAGCTGAGGGCAGCAG
TAAATCCAGGCCCGAATGGAACAGCAGCCCACTCGTCCTCCACAGAGCTCACAGCCACCACCACCTCCA
CCACCTATGCCATTAGAGCTCCAACGAAGCCTCCAGTTGGACCCAAAACCTTCTCCCTTGAAAGATAACC
CCTCACCTGAACCTCAGTTGGATGACATCAAAAGAGAGCTGAGGGCTGAAGTTGACATTATTGAACAAAT
GAGCAGCAGCAGTGGGAGCAGCTCTTCAGACTCTGAGAGCTCTTCGGGAAGTGATGACGATAGCTCCAGC
AGTGGAGGCGAGGACAATGGCCAGCCTCTCTCCGCAGCCTTCACACCAGCAGCCCTACAACAGTAGGC
CTGCCGTTGCCAATGGAACAGCCGCCACAAGGAAGCAACAGCTCATGAACGCCCTCAGAAATGACTT
GCAGTTGAGTGAGTCTGGCAGTGACAGTGATGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207169 protein sequence Red=Cloning site Green=Tags(s) MNGTANPLLDREEHCLRLGESFEKRPRASFHTIRYDFKPAIDTSCEGELQVGKGDEVITLPHIPGSTP PMTVFKGNKRPYQKDCVLIINHDTGEFVLEKLSSSIQVKKTRAEGSSKIQARMEQQPTRPPQTSQPPPPP PPMPFRAPTKPPVGPKTSPKDNPSPEPQLDDIKRELRAEVDIEQMSSSSGSSSSDSESSSGSDDDDSS SGGEDNGPASPPQPSHQPPYNSRPAVANGTSRPQGSNQLMNALRNDLQLSESGSDSD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_033083
ORF Size:	804 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_033083.6</u>
RefSeq Size:	4503 bp
RefSeq ORF:	807 bp
Locus ID:	85403
UniProt ID:	<u>Q96JC9</u>
Cytogenetics:	3p25.1
MW:	29 kDa